
Gebrauchsanweisung | Instructions for use |
Mode d'emploi | Modo de empleo | Modalità d'uso

CE 0483

remanium® 

Edelmetallfreie Legierungen für die Kronen- und Brückentechnik
Non-precious metal alloys for the crown and bridge technique
Alliages exempts de métaux précieux pour couronnes et bridges
Aleaciones exentas de metales preciosos para la técnica de coronas y puentes
Leghe non preziose per ponti e corone

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Dear Customer

Thank you for having chosen a quality product from Dentaureum.

In order to use this product at its best for you and your patients, it is important to study and follow these directions for use carefully.

The written instructions cannot cover all eventualities during operation. For this reason our Hotline is available to answer any other questions and ideas that may arise.

Due to a constant development we recommend, even when you use the same products frequently, to study the enclosed latest issue of the instructions for use or refer to our website at www.dentaureum.com

Manufacturer

Dentaureum GmbH & Co. KG | Turnstr. 31 | 75228 Ispringen | Germany

With reference to our quality standards

Dentaureum ensures a faultless quality of the products manufactured by us. These recommendations are based upon our own experiences. The user himself is responsible for the processing of the products. Responsibility for failures cannot be taken, as we have no influence on the processing on site.

Contents

1. General instructions for use	14
2. Waxing-up.....	14
3. Spruing	15
4. Calculation of metal quantity.....	15
5. Investing.....	15
6. Wax burn-out and preheating.....	16
7. Melting and casting	16
7.1. Melting crucible	16
7.2. Reuse of casting cones.....	16
7.3. Cast-on	16
7.4. Casting equipment.....	16
8. Induction melting with vacuum pressure casting machines.....	17
8.1. Premelting.....	17
8.2. Melting	17
9. Induction melting with centrifugal casting machines.....	17
10. Melting with torch	17-18
10.1. Observing the melt during torch casting.....	18
11. Overheating	18
12. Cooling of the casting ring	18
13. Devesting and sandblasting	18
14. Finishing	18-19
15. Processing and polishing steps.....	19
16. Preparation of surfaces for ceramic coating	19
17. Oxide firing	20
18. Ceramic veneering	20
19. Soldering	20
19.1. Soldering after firing	21
20. Laser welding.....	21
21. Cleaning.....	21
22. Contraindications and side effects.....	21

1. General instructions for use

These are the instructions for use when working with non-precious remanium® crown and bridge alloys.


Any deviations to the norm or specialized additional information about our non-precious remanium® crown and bridge alloys can be found in the specific alloy supplement supplied in each individual alloy package or in the internet under www.dentaurum.de.

Our dental technical team is available to answer any questions that may arise. Dental Technical Hotline Tel No: +49 72 31/803-410.

Information on contraindications and side effects can be found at the end of these instructions for use.

2. Waxing-up

The wax pattern is fabricated in the usual manner according to standard dental technical rules. Please only use waxes or modelling materials which burn out fully without leaving residues e.g. StarWax from Dentaurum.

Before waxing-up, the dies should be coated with a die spacer. The shape of the metal copings should correspond in a reduced scale to the final shape of the finished crowns. Missing portions of teeth should be compensated like fig. ① to ⑧.

In order to ensure that the metal fills the investment mould completely, the copings must be waxed to a minimum thickness of 0.4 mm.

A uniform thickness of ceramic ensures stress free adhesion to metal, see fig. ① to ⑧.


Apply ceramic to the framework in an equal thickness. The ceramic layer should cover the incisal edge of the crown in order to prevent chipping.

Areas to be soldered after casting (or firing), should be waxed up spaciouly, see fig. ⑤, ⑦ and ⑧. Rounding the metal/ceramic interfaces ① to ⑦ will provide neat junctures and aesthetic shading.

Pontics must be designed with either a collar or at least interdental reinforcement (Thermal conduction).


The outstanding properties of non-precious remanium® crown and bridge alloys allow very thin connections to be used between the crowns and pontics. Thus the teeth can be shaped more individually and a "block appearance" can be avoided.

3. Spruing

For bridges: Wax pattern with bar sprues

Casting button sprue	ø 3.5–4 mm
Runner bar	ø 4–5 mm
Connection to the crowns	ø 2.5–3 mm length 2.5–3 mm

Use the direct method for single crowns and small bridges

Sprue	ø 3 mm length 15–20 mm
Voluminous crowns	ø 3.5–4 mm length 15–20 mm

4. Calculation of metal quantity

General Rule: Wax weight of pattern including casting sprues multiplied by (8.2–8.6 g/cm³) plus 1 to 2 cylinder (approx. 10 g) for the casting button.

5. Investing

Special investment material for non-precious alloys e.g. rema® CC, Castorit® all speed, Castorit®-super C or Trivest (Dentaurum) with a high expansion.

A casting ring system using corresponding ceramic ring liners is recommended in order to have exact control over the final fit. "Kera-Vlies®" in 1 or 2 mm thickness is especially recommended.

6. Wax burn-out and preheating

For explicit information concerning conventional and speed preheating of the investment material please see the relevant manufacturer's instructions for use.

In general the final temperature for all remanium® alloys should be 950 °C / 1742 °F.

This standard temperature however, may vary according to the type of casting system used. The final temperature should be held for between 30 and 60 minutes according to the size of the casting ring.

7. Melting and casting

7.1. Melting crucible

Use a ceramic based crucible (i.e. magnesium oxide, silicon oxide or aluminium oxide).

Do not use carbon crucibles.

Preheat empty melting crucible in the furnace. Use one crucible per alloy and clean after each casting.

Do not use special flux powder!

For larger quantities of metal, a preliminary melting phase (casting delay phase) is recommended.

7.2. Reuse of casting cones

When the alloy is melted, some of the elements which later form the adhesive oxide are lost. If the alloy (e.g. in the form of a casting button) is re-melted several times, these elements may be diminished until they are completely destroyed. If the metal framework is later to be coated with ceramic, then only new material should be used for casting.

7.3. Cast-on

Metals used for cast-on and high platinum content attachments cannot be used in combination with non-precious remanium® crown and bridge alloys.

7.4. Casting equipment

Melting and casting techniques suitable for casting non-precious remanium® crown and bridge alloys are those which are capable of reaching the required casting temperatures. These include induction casting machines with vacuum-pressure or centrifuge, arc melting with vacuum pressure casting, or torch melting with Propane/Oxygen or Acetylene/Oxygen or other gas mixtures strong enough to reach the casting temperature.

8. Induction melting with vacuum pressure casting machines

Set vacuum to 250–450 mbar (recommended pressure: 450 mbar).

8.1. Premelting

Premelt the ingots until they are red hot (no visible rounding at the corners of the ingots).

8.2. Melting

Melt until the oxide layer of the molten metal breaks open, then cast immediately.

With exception of remanium® LFC: after the molten shadow disappears, continue heating for 2–4 seconds before the oxide layer cracks up, and then cast (see alloy info supplement).

Observe the melting procedure and initiate the casting procedure manually. Oxides formed during melting (cf. melting illustrations in the operating instructions for the pressure casting machine) remain as slag in the crucible after casting and can be easily removed. If the correct melting and casting technique is used, there is only a minimum residue in the crucible and the molten metal is not overheated.

Automatic casting with a pyrometer control is unsuitable, as the pyrometer cannot always accurately determine the right time for casting.

9. Induction melting with centrifugal casting machines

Wait until the last casting cylinder has collapsed, the oxide layer of the molten metal has broken and then cast immediately.


With exception of remanium® LFC: after the melt shadow disappears, cast immediately (see alloy info supplement).

We recommend premelting the metal if a large amount is to be cast.

10. Melting with torch

When melting, the remanium® crown and bridge alloys show a clean melt without sparks.

Caution! Use torch with shower head nozzle!

a) Propane/Oxygen Melting

Select maximum burner setting as stated by manufacturer.
Oxygen 2–3 bar.

b) Acetylene/Oxygen Melting

Follow instructions of burner manufacturer.
Acetylene: 0.7 bar, oxygen 3 bar.

The blue core of the flame at the burner head should be 4–5 mm in length. Pre-heat the empty crucible. Melt the metal in the crucible with a loose circular movement of the flame. The distance at which the torch nozzle should be held depends upon the type of nozzle and the gas mixture used.

10.1. Observing the melt during torch casting

When melting with a torch, an oxide layer forms when the last metal piece has collapsed. With a loose circular motion of the flame, continue the melting process until the metal under the oxide skin is seen to move under the pressure of the flame. Begin the casting process **before** the oxide layer cracks open.

11. Overheating

As with all melting and casting techniques, **overheating the metal must be avoided**. Overheating results in blowholes, micro porosities and grain coarsening. These casting imperfections are frequently reasons for cracking of a restoration or cracks in the ceramics.

12. Cooling of the casting ring

Remove the casting ring carefully and allow to bench cool.

Caution! Danger of fire! Always place hot casting rings onto heat resistant surfaces.

Ensure a safe distance from flammable objects.

13. Devesting and sandblasting

Avoid dust inhalation!

Moisten the casting ring before devesting in order to reduce the amount of dust produced!

Do not strike the casting button with a hammer!

Sandblast using pure aluminium oxide, grain size approx. 125 µm.

14. Finishing

General information:

Caution: Always use suction units when grinding, cutting and polishing!

The higher physical characteristics of the non-precious alloys require grinding and polishing tools different to those normally used for processing precious metal alloys.

The grinding tools selected should only be used for one alloy.

Observe the minimal coping thickness of 0.2–0.3 mm.

Due to the high physical properties of non-precious alloys, both the wall thickness of the crown and the thickness of the bridge connecting points can be reduced by approx. 30 % in comparison to average precious metal alloys.

Always grind in one direction without applying too much pressure. Avoid overlapping of the metal. Remove irregularities such as porosities and entrapments.

15. Processing and polishing steps

Steps	Tools	REF	Remarks
Removal of casting sprues	Supercut STM separating discs	130-111-00 130-113-00	High cutting power and long life span. Universal in application.
Processing, rough	Grinding wheel	131-322-50	Grinding of sprue attachments. Rough processing, e. g. of interdental spaces (surfaces intended for ceramic coating must be processed with tungsten carbide tools afterwards).
	Supercut STM separating discs	130-112-00	
Processing, fine	Tungsten carbide burs	123-582-00 123-584-00 123-585-00 123-601-00	General processing with "rough" tungsten carbide burs. Fine grinding can be done with fine tungsten carbide burs.
Fine grinding	Aloxin stones	135-852-00 135-853-00	Pure aluminium oxide abrasives for smooth even transitions, especially metal occlusal surfaces.
Buffing	Silichrom polishers	138-645-00 138-640-00	Fast effective material removal (rough buffing).
	Grey rubber polishers	138-102-00 138-302-00	Medium abrasive effect.
	Green rubber polishers	138-101-00 138-301-00	Fine abrasive effect.
Polishing	Polishing brushes	141-800-00	Universal polish with hand tool.
	Polishing paste Tiger brillant	190-350-00	Fast acting high shine polishing paste.
	Universal Finishing Paste Tiger Starshine	190-301-00	The ideal paste for final finishing.
Burnishing of crown's inner surfaces	Al ₂ O ₃ – 50 µm, high lustre shot blasting beads	128-017-00 128-211-00	For use with pen-type blasting tool. Important! Cover ceramic edges with wax.

16. Preparation of surfaces for ceramic veneering

Grind the metal surfaces with tungsten carbide tools to ensure a gradual transition. Blast the surfaces with a fine "pen-type" shot blasting tool using pure aluminium oxide of medium grain (125 µm) at a low pressure of 2–3 bar. Clean with ultrasonic cleaner in distilled water.

17. Oxide firing

Oxidation firing is not required with remanium® bonding alloys. It is advisable to fire an oxide bake in order to visually check the framework condition (5 minutes without vacuum at the firing temperature of the opaque, unless otherwise indicated by the ceramic manufacturer). Following the bake, the oxide should be removed by blasting with single use aluminium oxide abrasive, grain size 125 µm and low pressure of 2–3 bar. Afterwards the surface must be cleaned again.

18. Ceramic veneering

remanium® ceramic bonding alloys can be veneered using conventional ceramics such as ceraMotion® Me or ceraMotion® LFC (Dentaurum) as long as the CTE is suitable with the alloy (see alloy info supplement).

The coefficient of thermal expansion (CTE) of many ceramic varieties increases, depending on the amount of times the ceramic is fired and the length of time within the firing furnace.

The longer and the more often the ceramic is fired, the higher the CTE. The metal CTE remains constant; the adhesion takes place by applying cooling times after firing and standby temperatures.

The opaque material can be applied to the sandblasted and cleaned metal framework according to the ceramic manufacturer's recommendation.

It is important to dry the opaque thoroughly. Always ensure the ceramic furnace is calibrated!

The ceramic is then applied according to the instructions for use!

If no other information is given by the ceramic manufacturer, the object is cooled as recommended in the alloy info supplement.

Important: After each firing, brush the framework under running water and allow it to dry.

19. Soldering

Try to avoid soldering if possible in order to prevent a material mixture. If however, it is necessary to solder, please ensure a suitable solder with the correct composition and melting temperature which is appropriate for the alloy used.

Roughen the joint areas with fine corundum stones.

Connect the parts using wax or acrylic, remove from the dies and fix in investment material for soldering. Boil out. Apply a coat of flux to the solder joint areas. Allow the model to dry, then preheat at 600 °C/1112 °F for 10 min. Apply a further coating of flux and heat the casting until it is red hot. Cut the solder to the appropriate length, dip the pieces into the flux and place them into the soldering gap. The flux on the pieces of solder and soldering gap must be heated with a hot flame until it covers the joint area completely. Only then bring the solder to flow through the hot flame.

Caution: Soldered units should cool slowly after the ceramic has been fired.

19.1. Soldering after firing

It is not recommended to solder after having applied the ceramic. The reason for this is that the corrosive resistance is reduced and the low diffusion ability of precious metal solders to non-precious alloys. Recommended alternatives are joining techniques such as laser welding, WIG welding and adhesive joining.

20. Laser welding

The laser welding technique produces a solder-free, mechanically strong and corrosion resistant joint.

It is important to observe the framework design, the surface structures, the welding sequence and the welding parameter of each individual laser machine. Suitable welding material is welding wire in the same metal as the remanium® crown and bridge alloy.

21. Cleaning

Crowns and bridgework made in remanium® crown and bridge alloys can be cleaned, after finishing, with ultrasonic solution P.

22. Contraindications and side effects

Signs of intolerance to non-precious remanium® crown and bridge alloys after are rare correct processing.

If the patient has a proven allergy against any component within the alloy, for safety reasons this alloy must not be used.

If various alloys have been used, it is possible that galvanic effects may occur.

Some patient have shown traces of local irritation caused by electrochemical reactions.

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