Orthocryl®

CE 0483

| | |
|----------------------|-------|
| Gebrauchsanweisung | S. 2 |
| Instructions for use | P. 10 |
| Mode d'emploi | P. 18 |
| Modo de empleo | P. 26 |
| Modalità d'uso | P. 34 |

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Dear Customer

Thank you for choosing this high-quality product from Dentaaurum. In order to achieve the best results with this product, it is important to carefully study and follow these instructions for use.

The written instructions cannot cover all possible situations that can occur during use. For this reason, our hotline is available to answer any other questions or concerns that may arise.

Due to constant developments, we recommend, even when using the same products regularly, that you study the enclosed current instructions for use or refer to our website at www.dentaaurum.com.

1. Manufacturer

Dentaaurum GmbH & Co. KG
Turnstr. 31 | 75228 Ispringen | Germany

2. General product description

Areas of application

Orthocryl®/Orthocryl® EQ are PMMA-based, cold-curing polymers for the fabrication of fixed and removable orthodontic appliances, bite-guards, anti-snoring devices etc.

All Orthocryl® powders (polymers), liquids (monomers) and colour concentrates are Class IIa medical devices and are approved materials for use in fixed orthodontic techniques.

Accessories such as Orthocryl® Disco Glitter and novel designs are intended solely for use in removable appliances.

Orthocryl® achieves optimal stability with the spray-on technique, but it is also suitable for the doughing technique.

Orthocryl® EQ achieves optimal flowing characteristics in the doughing technique, but it is also suitable for the spray-on technique.

3. Application

Firstly the models must be „**de-aerated**“. Otherwise the air bubbles emerging from dried plaster models could cause the acrylic to lift from the model, or penetrate into the resin! Therefore the models should be immersed in warm water – depending on plaster quality – for approximately 10-25 min, until no more air bubbles emerge from the plaster.

In order to avoid direct skin contact whilst processing Orthocryl® components, it is mandatory to wear facial masks, protective goggles and special protective gloves such as Dermatril® disposable protective gloves made from special-nitrile (REF 162-907-00,162-908-00 and 162-909-00). For further information, please see the relevant safety data sheets.

In order to achieve smooth, silky-matt, polished basal appliance surfaces, without discolouration or plaster inclusions, we recommend coating the model after soaking with a special Orthocryl® separating medium (REF 162-800-00).

3.1. Applying the Orthocryl® separating medium (REF 162-800-00)

- Carefully blow the soaked model.
- Pour the separating medium directly from the bottle (REF 162-100-00) onto the moist model and spread evenly over the surface using a paintbrush.
- Stand the model onto its back so that any excess fluid may drain away.

- After approx. 30 seconds, carefully blow away separation fluid residues or droplets with the compressed air and remove from the screw and wire retention if necessary.

Begin with the acrylic build-up using the spray-on (salt and pepper) technique (3.2) or the doughing technique (3.3).

3.2. Spray-on (salt and pepper) technique

First, apply the powder (polymer) then the liquid (monomer) directly onto the prepared model, alternating with gentle pendulum movements.

Only use as much liquid as the powder can absorb. The resin should not slump!

In order to avoid air bubbles around the screw, use more liquid in this area at the beginning in order to allow the acrylic to flow under the screw without hindrance.

Finally apply enough powder so that a dry powder layer remains. This should assure an optimal fit with minimal shrinkage.

3.3. Doughing technique

Mix the Orthocryl®/Orthocryl® EQ in a ratio of 2.5 : 1

(Volume parts – powder : liquid). Allow time to swell, with Orthocryl® approx. 5-7 min and with Orthocryl® EQ approx. 3-5 min. Times may be shorter at higher room temperatures. Then apply the acrylic directly to the model.

3.4. Application of Orthocryl® color concentrates

The Orthocryl® liquids can be individually dyed with Orthocryl® color concentrates. The recommended mixing ratio is between 1 : 10 to 1 : 50 (Color concentrate : Orthocryl® clear liquid). The dosage pipette (REF 162-101-00) guarantees accurate dosage of the color concentrate.

For the standard colours we recommend:

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3.5. Application of Orthocryl® Disco Glitter

Orthocryl® Disco Glitter can be mixed together with the Orthocryl® powder as much or as little as desired. Neither the application characteristics nor the properties of Orthocryl® are influenced by this. We recommend mixing only a little glitter, for example one heaped knife edge or Le Crons. Mixing ratio: 0.5 : 100 (ratio of glitter : powder).

3.6. Application of novel designs

The novel designs – for jazzing up appliances – look best in very transparent acrylics. As the novel designs do not chemically bond with the acrylic, these must be completely covered with acrylic.

3.7. Processing Orthocryl® black & white

Please note, that the Orthocryl® „black & white“ acrylic has a very high degree of opacity. Orthocryl® Disco Glitter and the picture designs lose their effect due to the high opacity of the acrylic.

When finishing, it is important to avoid grinding the screws and wires. The opacity of Orthocryl® „black & white“ can be reduced by adding clear Orthocryl® acrylic.

A mixing ratio of 50% clear polymer and 50% “black” or “white” polymer reduces the opacity, so that wires, screws etc. become visible in Orthocryl®, while the acrylic still exhibits a white or black appearance.

You can also choose other mixing ratios to mix the polymers and to obtain more or less opacity.

4. Curing/Polymerization time

Pressure pot: 2.2 bar/32 psi

Temperature: 40-46°C/104-114.8° F

Time: No less than 15 minutes (Orthocryl® EQ)!

No less than 20 minutes (Orthocryl®)!

Polymerisation in a pressure vessel should not be interrupted during the first 8 minutes (Orthocryl® EQ) or 10 minutes (Orthocryl®) in order to prevent micro porosity from occurring.

Please observe the directions for use and safety instructions of the pressure pot!

5. Finishing and polishing

Orthocryl® acrylics can be finished and polished using conventional acrylic burs and polishers. Please observe the stipulated protective equipment which must be worn (facial protection/protective eyewear) during processing.

6. Reducing the remaining monomer content and the biocompatibility of Orthocryl® acrylics

In order to reduce the remaining monomer content, the finished appliance – if possible – should be immersed in water for three days, before transferring to the patient. This immersion in water causes the monomer residues to reduce to values similar to heat-curing resins. If this is not possible, then it is recommended that the patient places the appliance in water, during the time it is not worn. Orthocryl® acrylics have been scientifically tested.

Information on this, a comprehensive instruction manual with many practical tips and the entire product range of Orthocryl® is available for your use.

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7. Instructions for cleaning and disinfection

Orthocryl® acrylics should be cleaned with tooth paste and brush. Orthocryl® acrylics can be disinfected with all commercially available disinfectants for impression materials (e.g. Impresept®/3M Espe, MD520/Dürr Dental, Dentavon®/Schülke+Mayr or Mucalgin®/Merz Hygiene).

Please observe the instructions by the respective manufacturer.

8. Information on the constituents

Monomer = Methyl methacrylate „MMA“ (liquid)

Polymer = Polymethyl methacrylate „PMMA“ (powder)

Please observe the safety data sheet which every new customer automatically receives on first time use and which is available upon request!

9. Warning and safety instructions

Danger contains methyl methacrylate and ethylene dimethacrylate
Causes skin irritation.

May cause an allergic skin reaction.

Causes serious eye irritation.

May cause respiratory irritation.

Keep away from heat/sparks/open flames/hot surfaces. - No smoking.

Wear protective gloves/protective clothing/eye protection/face protection.

Please observe the safety data sheets which every new customer automatically receives on first time use and which is available upon request!

If you are allergic to one or more constituents, please refrain from using this product!

10. Storage and service-life

Powder (Polymer): 10 years in a closed container

Liquid (Monomer): 36 months from the date of production, in a closed bottle

- Protect from contamination and sunlight as they may cause premature polymerisation.

11. Disposal

Please observe the requirements stated in the relevant safety data sheets for the disposal of non-polymerised Orthocryl® components.

12. Tips for cleaning and caring for Orthocryl® appliances

Appliances should be cleaned regularly and stored in water when not in use to prevent problems with hygiene. This procedure prevents crack corrosion or brown discoloration of the wire components.

13. Delivery

The complete range of Orthocryl® products can be found in the orthodontics or prosthetics catalogue.

14. With reference to our quality standards

Dentaurum ensures a faultless quality of the products that we manufacture. These recommendations are based upon our own

experiences. The user himself is responsible for the processing of the products. Responsibility for failures cannot be taken, as we have no influence on the processing on site.

Please refer to the label. Additional information can be found at www.dentaurum.com (Explanation of the symbols REF 989-313-00).

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