DE

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2. Brief description

tiologic® implant types are designed for insertion in the endosteal region of the maxilla or mandible. Depending on the indication, corresponding transgingival abutments are secured on the implants and fitted with a prosthetic superstructure.

The tiologic® implant system contains specially coordinated instruments, abutments and accessories for placement of the implants and fabrication of the prosthetic restoration. Only original components of the tiologic® implant system should be combined in accordance with the instructions for use / user manuals.

3. Indications

tiologic® implant types can be used in both the mandible and maxilla for surgical immediate implantation, delayed immediate implantation and delayed implantation using either the one-stage or two-stage technique. Indications for implant placement are small and large bounded saddles (one-tooth restorations, increasing the number of abutments) in the maxilla and mandible, a shortened dentition and an edentulous jaw. Before the decision, the risks and advantages / disadvantages of implantation as well as alternative methods of treatment should be considered.

With all implantological cases the implant diameter and length of the tiologic® implant types should be in proportion to the prosthetic restoration. Implants with a minimum diameter of 4.2 mm should always be used for restorations that subject the implant and superstructure to high mechanical loading, if this is practical with the particular oral situation.

We offer tiologic® S ø 3.3 mm implants for patients with narrow alveolar ridges. Due to the small diameter and low load capacity (compared to the tiologic® M ø 4.2 mm implants), these implants have a limited range of indications. In fully edentulous cases, four or more tiologic® S ø 3.3 mm implants must be inserted with a splinted bar restoration without extension. In partially edentulous cases, implant supported restorations must be combined with tiologic® ø 4.2 mm, ø 4.8 mm or ø 5.5 mm implants and a splinted fixed prosthetic restoration. In single restorations, tiologic® S ø 3.3 mm implants should only be used for the lower incisors or the upper lateral incisors and only with a minimum 11.0 mm implant length. Single restorations on tiologic® ø 3.7 mm, ø 4.2 mm, ø 4.8 mm or ø 5.5 mm implant types require a minimum 9.0 mm implant length. Care should be taken to avoid an excessive mechanical loading when using ball head abutments together with ø 3.3 mm implants.

4. Contraindications

Implants with a diameter of 3.3 mm are not suitable for single-tooth restorations of the central incisor in the maxilla or the canines, premolars or molars in the maxilla or the mandible. It is not permitted to use telescope crown constructions on these implants.

The use of LOCATOR® abutments for non-parallel abutments of 10° or more per implant is contraindicative.

General contraindications:

General contraindications for dental surgery procedures apply. These include: reduced immunodeficiency / steroid treatment / blood coagulation disorders / uncontrolled endocrine diseases / rheumatic disorders / bone system diseases / cirrhosis of the liver / drug, alcohol or tobacco abuse / depression, psychopathic disorders / poor patient compliance / chronic inflammatory underlying diseases.

Local / personal contraindications:

Osteomyelitis / radiotherapy in the head region / recurring mucosal diseases / temporomandibular joint dysfunctions / parafunctions / lack of vertical or horizontal bone availability, jaw defects, inadequate bone quality / poor oral hygiene.

5. Further information

Though placement of dental implants has a high rate of success and implants have a long durability, successful treatment cannot be guaranteed. The operator should note and document any problem cases and inform the manufacturer.

An inadequate number of implants, implants with insufficient length or diameter, unfavorable positioning of the implants or a statically poor prosthetic restoration can cause fatigue fracture in implants, abutments and prosthetic screws under bio-mechanical loading. Placement of the implants and fabrication of the prosthetic restoration should take into account the particular oral situation to avoid overloading the components.

Using tiologic® implant system components in combination contrary to the instructions for use / manuals and which are not described as being specifically combinable, can also cause mechanical failure, damage to the tissue or unsatisfactory aesthetic results.

tiologic® implant types are not known to have any side effects or produce any interaction. It cannot, however, be ruled out that in rare cases allergies to components used in the materials of the tiologic® implant system may occur or that there may be electrochemically-induced discomfort.

6. Packaging and sterility

All tiologic® implant types are supplied singly with the appropriate closure screws in gamma-sterilized double packaging(foil and blister packaging) that protects the content against contamination. Outer packaging provides added protection to the double packaging. If the double packaging is visibly damaged, the product must not be used. Outer packaging provides added protection to the double packaging. The tiologic® implants are single-use products.

All other components of the tiologic® implant system are supplied non-sterile.

The surgical components should be thoroughly disinfected, cleaned and sterilized before use. This applies in particular to components that can be disassembled, e.g. the torque wrench.

7. Application / Availability / Precautions / Documentation

The tiologic® product range is supplied exclusively to doctors, dentists and dental technicians. It should only be used by doctors, dentists or dental technicians who are familiar with dental implantological procedures, including diagnosis, preoperative planning, surgical technique and prosthetic treatment.

Before use, operators should ensure that they have carefully read and understood the full tiologic® instructions for use / manuals. This quick reference guide does not replace the instructions for use manuals for the tiologic® implant system. As the instructions and manuals cannot provide all information for immediate use, we strongly recommend that, before using the system, operators attend a tiologic® system training course offered by Dentaurum Implants to learn the correct techniques.

- Refer to the Product Catalog and the Surgery Manuals for information on precautions and the selection of components for the clinical procedure.

- Refer to the Product Catalog and the Prosthetic Manuals for information on precautions and the selection of components for the prosthetic procedure.

Before using this product, the operator must give the patient a thorough examination and a detailed explanation of the procedure. Dentaurum Implants recommends full clinical, radiological, photographic and statistical documentation.

The tiologic® implant system components can be documented, e. g. in the patient file or PatientPass, using the additional labels.

The operator should ensure the products cannot be aspirated during intraoral use.

Not all components are available in every country.

8. Quality / Warranty and Liability / Development

Development, clinical testing, production and quality control of the tiologic® product range are completed in accordance with the Medical Device Directive 93/42/EEC.

Sections 9 and 10 of our General Terms of Delivery and Payment apply with regard to warranty or liability – unless otherwise provided in this quick reference guide.

Warranty and liability are rendered void in particular if the products are not used by the operator or a third party in accordance with the instructions for use; this also applies if the tiologic® product range is used in combination with products from other manufacturers, which have not been specifically recommended for use by Dentaurum Implants.

Dentaurum Implants has no control over processing and application of the product. These are the sole responsibility of the user.

Technical advice (oral and written) is based on the scientific and technical knowledge available when the product is put on the market. It does not release the user from the responsibility of personally checking the suitability of the products for the intended indication and application. Advice is only given as a non-binding recommendation, which cannot be assumed to provide any form of assurance or guarantee.

All products are subject to continuous development based on current scientific knowledge and we reserve the right to make changes in the construction, design or material of the products.

9. Explanation of symbols used on the label

Please refer to the label. Additional information can be found at www.dentaurum-implants.com (Explanation of the symbols used on the label REF 989-313-00).

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