DE

Gebrauchsanweisung
RESOLON® (PA 6-6/6)
steriles, nicht-resorbierbares
Nahtmaterial

EN

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RESOLON®

EN

Instructions for use, please read carefully!

The information given in this package leaflet is updated regularly. Please read instructions carefully prior to use.

Polyamide 6-6/6

sterile, nonabsorbable
suture material

DESCRIPTION

RESOLON® is a monofilament, non-absorbable, surgical suture material made from polyamide 6-6/6, extruded from a copolymer of polyamide 6 and polyamide 6/6. RESOLON® is available dyed blue (copper phthalocyanine blue, Color Index 74160).

RESOLON® is available in various thicknesses and lengths and also as a needle-thread combination. RESOLON® satisfies the requirements of the European Pharmacopoeia for sterile nonabsorbable suture material.

INDICATIONS

RESOLON® is generally used for sutures and ligatures in soft tissues including use in neurosurgery and ophthalmic and cardiovascular surgery.

USE

The suture material should be chosen in accordance with the condition of the patient, the experience of the surgeon, the surgical technique, and the size of the wound.

ACTION

RESOLON® causes only a slight inflammatory tissue reaction initially and is then slowly and progressively encapsulated by connective tissue. Although RESOLON® is not absorbed, gradual hydrolysis of the polyamide in vivo can eventually lead to a progressive loss of tensile strength.

CONTRAINDICATIONS

As RESOLON® eventually loses some of its tensile strength, it is not suitable for indications in which lasting stability of the suture material is required.

WARNINGS / PRECAUTIONS / INTERACTIONS

As with all suture material, care must be taken to ensure that the thread is not damaged during handling. In particular, it must not be kinked or crushed by surgical instruments such as needle holders.

When tightening the suture always pull on the thread between the needle and the puncture channel. Do not pull the thread too strongly or over sharp objects. When stretching the thread avoid friction with the rubber glove, as this can damage the thread.

In order not to damage the needle, always grasp it 1/3 – 1/2 of the distance from the end to the point. Do not bend the needle, as this leads to a loss of stability. Because of the risk of infection, the user should take particular care not to incur stab wounds when using surgical needles. Used needles must be disposed of correctly.

The risk of wound dehiscence varies with the location of the wound and the suture material used, therefore the user should be familiar with the surgical techniques in which RESOLON® is to be used. As with all foreign bodies, prolonged contact between RESOLON® and saline solutions can lead to the formation of calculi (urinary and biliary tracts).

RESOLON® is used where indicated in accordance with standard surgical suturing and knotting techniques and the experience of the user.

ADVERSE REACTIONS

Inflammatory tissue reactions and transient irritation in the area of the wound may occur initially. The foreignbody effect can exacerbate pre-existing tissue inflammation initially.

STERILITY

RESOLON® is sterilised by ethylene oxide or gamma radiation. Do not use opened or damaged packs! Do not resterilise!

STORAGE CONDITIONS

Store at not more than 25°C protected from moisture and direct heat. Do not use after the expiry date!

DESCRIPTION OF THE SYMBOLS USED ON THE PACKAGE

REF = Reference number

LOT = Batch number

 = Use by year - month

 = Consult instructions for use

 = Do not reuse

 = Do not resterilize

 = Do not use if package is damaged

STERILE EO = Sterilised using ethylene oxide

STERILE R = Sterilised using irradiation

 ^{25°C} = Upper limit of temperature

CE ₁₂₇₅ = CE marking and identification number of the notified body. Product conforms to the essential requirements of the Council Directive 93/42/EEC concerning medical devices.

 = dyed, monofilament, nonabsorbable

PA = Polyamide

 = HIBC code

 = loop suture

 = ligature pack:
sterile individual suture on a small roll

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RESORBA Wundversorgung GmbH & Co. KG